

Essential Fitting Guide

Lenses>SpecialtyGP>Multifocal> ESSentialGP>Fitting Guide

Parameter Selection:

1. Initial Base Curve/Diameter Selection

Select initial base curve according to the base curve selection chart.

Select diameter according to diameter selection chart.

Base Curve Selection Chart

Select base curve according to corneal cylinder

Determine Flat K	0.00D to 0.62D	0.75D and UP
40.00 to 40.37	8.30	8.20
40.50 to 40.87	8.20	8.10
41.00 to 41.37	8.10	8.00
41.50 to 41.87	8.00	7.90
42.00 to 42.37	7.90	7.80
42.50 to 43.00	7.80	7.70
43.12 to 43.62	7.70	7.60

43.75 to 44.25	7.60	7.50
44.37 to 44.87	7.50	7.40
45.00 to 45.50	7.40	7.30
45.62 to 46.12	7.30	7.20
46.25 to 46.75	7.20	7.10
46.87 to 47.37	7.10	7.00

Diameter selection chart

Select diameter according to corneal cylinder

Base curve	Myopes	Hyperopes
$\leq 7.30\text{mm}$	9.0mm	9.2mm
7.40 to 7.90mm	9.2mm	9.5mm
$\geq 8.00\text{mm}$	9.5mm	9.5mm

2. Add Series Selection

Select add series according to patient's add

Patient's Add	Add series
+0.50 to +1.50	Series 1
+1.75 to +2.25	Series 2
+2.50 and up	Series 3

3. Lens Position and Movement

Evaluate lens position and movement; the ideal fit will be superior central (upper lid attachment) with a fluorescein pattern that demonstrates alignment along the flattest corneal meridian. Make base curve and diameter changes accordingly (see Troubleshooting guide).

4. Lens Power

Perform your over-refraction with loose trial lenses to determine the final distance Rx. Expect final Rx to be -0.50 D more than the existing contact lens Rx. Place the over-refraction in a trial frame and evaluate the transition from distance to near vision. If the over-refraction leads to acceptable distance but unacceptable near vision, reassess your base curve and/or add selection.

I Important Note: *In order to maximize the ADD available in each series, the lens needs to translate upward along the vertical corneal meridian as the patient looks from distance to reading tasks. An upper lid attachment will facilitate the upward transition of the lens.*

Troubleshooting Guide

Clinical Observations	Common Symptoms	Fitting Goal	Remedy
High-Riding Lens	Excellent near, Excessive minus Over-refraction to	Lower the lens to pick up more distance power in the center of the lens.	Steepen B.C. by 0.10mm and specify a thin edge design

achieve good

distance.

And/Or

Reduce lens diameter by 0.3mm, and specify a thin edge design

Last Alternative

Steepen B.C. by 0.10, reduce diameter by 0.3mm, specify a thin edge design and increase the lens center thickness by 0.05mm

Lens Rides Temporally

Near is usually excellent, distance requires additional -0.50D to -0.75D.

Get lens to center more by changing the dynamics between the upper lid and the lens.

Diameter \leq 9.3mm

Increase lens diameter by 0.50mm and specify a minus carrier.

Diameter > 9.3mm

Increase lens diameter by 0.3mm and specify a minus carrier

Low-riding lens

Good to excellent distance, near almost non-existent.

Need to achieve an upper lid attachment, making the greater plus power on the periphery of the lens available for near tasks.

Diameter \leq 9.3mm

Increase lens diameter by 0.5mm and specify a minus carrier

Diameter > 9.3mm

Increase lens diameter by

0.3mm and
specify a minus
carrier

Last Alternative

Flatten B.C. by
0.10mm (only if
the fit allows to do
so).

Lens that centers

Excellent
distance, reading
difficult

Get the lens to
translate upward
upon downward
gaze.

Increase lens
diameter by
0.3mm and
specify a minus
carrier.

Lens Rides Nasally

Distance is
difficult and near
vision is almost
non-existent.

Get the lens to
center more, by
changing the
dynamics
between the upper
lid and the lens.

**Diameter \leq
9.3mm**

Increase lens
diameter by
0.50mm and
specify a minus
carrier

Diameter $> 9.3\text{mm}$

Increase lens
diameter by
0.30mm and
specify a minus
carrier.