

# Tangent Streak Bifocal Fitting Guide

Lenses>Specialty GP>Bifocal>Buckley>Fitting Guide

Parameter Selection:

1. Select lens diameter based on flattest K and/or the horizontal visible iris diameter.

Visible Iris Diameter	Flat K Reading	Diameter
Less than 11.3	47.00D to 45.00D	9.1/8.8
11.3 to 11.8	44.75D to 43.00D	9.5/9.1
Greater than 11.8	42.75D to 40.00D	9.9/9.5

- For every .3mm increase in horizontal diameter, flatten base curve 0.25D.
- For every .3mm decrease in horizontal diameter, steepen base curve 0.25D.

2. Select base curve based on flattest K and amount of corneal astigmatism.

Flat K Reading	Amount of corneal astigmatism
----------------	-------------------------------

	0.00D	0.50D	1.00D	<1.00D
<b>47.00D</b>	47.00/7.18	47.50/7.10	48.00/7.03	See Below
<b>46.75D</b>	46.75/7.21	47.25/7.14	47.75/7.06	
<b>46.50D</b>	46.50/7.25	47.00/7.18	47.50/7.10	
<b>46.25D</b>	46.25/7.29	46.75/7.21	47.25/7.14	
<b>46.00D</b>	46.00/7.33	46.50/7.25	47.00/7.18	
<b>45.75D</b>	45.75/7.37	46.25/7.29	46.75/7.21	
<b>45.50D</b>	45.50/7.41	46.00/7.33	46.50/7.25	
<b>45.25D</b>	45.25/7.45	45.75/7.37	46.25/7.29	
<b>45.00D</b>	45.00/7.50	45.50/7.41	46.00/7.33	
<b>44.75D</b>	44.50/7.58	45.25/7.45	45.75/7.37	
<b>44.50D</b>	44.25/7.62	45.00/7.50	45.50/7.41	
<b>44.25D</b>	44.00/7.67	44.75/7.54	45.25/7.45	
<b>44.00D</b>	43.75/7.71	44.50/7.58	45.00/7.50	
<b>43.75D</b>	43.50/7.75	44.25/7.62	44.75/7.54	
<b>43.50D</b>	43.25/7.80	44.00/7.67	44.50/7.58	

<b>43.25D</b>	43.00/7.84	43.75/7.71	44.25/7.62
<b>43.00D</b>	42.75/7.89	43.50/7.75	44.00/7.67
<b>42.75D</b>	42.25/7.98	43.25/7.80	43.75/7.71
<b>42.50D</b>	42.00/8.03	43.00/7.84	43.50/7.75
<b>42.00D</b>	41.50/8.13	42.50/7.94	43.00/7.84
<b>41.75D</b>	41.75/8.08	42.25/7.98	42.75/7.89
<b>41.50D</b>	41.50/8.13	42.00/8.03	42.50/7.94
<b>41.25D</b>	41.25/8.18	41.75/8.08	42.25/7.98
<b>41.00D</b>	41.00/8.23	41.50/8.13	42.00/8.03
<b>40.75D</b>	40.75/8.28	41.25/8.18	41.75/8.13
<b>40.50D</b>	40.50/8.23	41.00/8.23	41.50/8.23
<b>40.25D</b>	40.25/8.38	40.75/8.28	41.25/8.18
<b>40.00D</b>	40.00/8.43	41.00/8.23	41.00/8.38

- Corneal toricity greater than a diopter, select a base curve steeper than flat K by a  $\frac{1}{4}$  of the amount of corneal toricity.

### 3. Power Calculation

### **Distance Power**

- The distance power is calculated by adding plus or minus to the spherical power of the manifest refraction depending on whether the base curve selection is flatter or steeper than the flattest K reading.
- For best results over refract with hand held lenses in normal lighting over the best fitting diagnostic lens.

### **Add Power**

- Use the spectacle add or for best results, use hand held trial lenses to over refract in adequate room illumination while patient is in the reading position

## **4. Seg Height**

- When visible white exists ( $< 1.00$  mm) between inferior limbus and lower eyelid, choose a 4.4 mm.
- When inferior limbus and lower lid are tangent, choose a 3.9 mm.
- When inferior limbus is above lower eyelid (not  $> 0.75$  mm above), choose 3.5 mm.

## **Lens Fitting**

### **A good fit exhibits the following:**

Intrapapillary fit

In distance gaze, the lens rests on the lower lid.

In near gaze, the lens has good upward translation.

The seg height is located at the lower pupil margin.

The edge lift is adequate.

The seg height is positioned at  $180^\circ$ .

The seg height drop is  $< 1.0$  second.

### **A flat fit is determined by the following objective findings:**

The lens decenters laterally following the blink.

There is excessive edge lift.

The lens rotates and rocks on the blink.

There is an absence of tears centrally.

The segment height angle fluctuates.

The lens drops extremely fast.

The lens decenters temporally.

**A steep fit is determined by the following objective findings:**

There is an excessive amount of pooling centrally.

The lens is bearing mid-peripherally.

There is no or very little edge lift.

The lens is decentering either inferiorly or nasally.

The segment drop is slow or not at all.

There is little or no translation.

Vertical movement is limited.

The segment has  $>15^\circ$  nasal rotation.

## Problem Solving Techniques

### Objective Findings

### Parameter Change

Superior positioning lens with normal edge lift

Steepen the base curve and/or an increase of prism by .50D.

Nasal or temporal dislocation

Larger diameter

The lens doesn't drop quickly enough.

The base curve may need to be adjusted either flatter or steeper (depending on the fluorescein pattern). Also consider increasing the prism

The patient holds their chin to the chest for distance vision.

Check the distance of the seg height in relationship to the pupil with the aid of the ophthalmoscope in primary gaze. It should be within the lower pupil measure. The seg height may need to be lowered.

The patient holds the reading material to their chest.

Raise the seg height.

The lens is rotating more than 15°, nasal.

Flatten the base curve by a .50D or/add .50D prism.

The lens is rotating more than 10°, temporal.

Steepen the base curve by .50D and/or add .50D prism.

The Lens pops out.

Steepen the base curve.

If the lens slips underneath the lower lid.

Increase the prism by .50D.

**Subjective Findings**

The patient observes superior flare.

**Parameter Change**

A larger diameter is needed.