



## FITTING GUIDE

### BASE CURVE SELECTION\*

CORNEAL CYLINDER	BASE CURVE
0.00 to 1.00	On Flat K
1.12 to 1.50	.25D steeper than Flat K
1.62 to 2.00	.50D steeper than Flat K
2.00 to 2.50	.75D steeper than Flat K
2.75 or more	discuss with consultant

\*This nomogram is based on using the recommended diameter.  
 Smaller Diameter - Fit .25D steeper  
 Larger Diameter - Fit .25D flatter

### DIAMETER SELECTION

BASE CURVE	DIAMETER RECOMMENDED
7.10 mm to 7.38 mm	9.0mm
7.42 mm to 8.08 mm	9.3mm
8.13 mm to 8.49 mm	9.6mm



Manufactured  
exclusively in

PARAGON HDS

## PARAMETER AVAILABILITY

32 Base Curves	7.10mm to 8.49mm
65 Powers	+6.00D to -10.00D
3 Diameters	9.0, 9.3, 9.6
Material	Paragon HDS *
Color	Blue

## THICKNESS CHART

Power (Diopters)	Thickness (mm)
Plano to -0.75	.16
-1.00 to -1.25	.15
-1.50 to -2.00	.14
-2.25 to -2.75	.13
-3.00 to -3.75	.12
-4.00 to -4.75	.11
-5.00 to -10.00	.10
+0.25 to +6.00 thicknesses available upon request.	



1.800.241.9312

**We Fit Your Practice™**