

MetroFocal Toric

Multifocal Soft Lens

LENS PARAMETERS*

Fitting Guide

| | |
|--------------|--------------------|
| Base Curve - | 8.6mm and 8.9mm |
| Power - | +10.00D to -10.00D |
| Cylinder - | -.75D to -6.00D |
| Axis - | Around the clock |
| Add Power - | up to +2.50D |
| Diameter - | 14.0mm |
| Color - | Light Blue |

* Custom parameters available at an additional charge.

Manufactured on advanced computerized lathing systems in Polymacon (38% water content material) for accuracy, optical clarity and guaranteed performance.

PATIENT SELECTION

Corneal cylinder from .75D to 6.00D
Needing up to +2.50 add

FITTING SYSTEM

(No trial lenses are required for fitting.)

With prescription in minus cylinder form, reduce sphere power by .25D.

Adjust cylinder power as follows:

Refractive Cylinder

-2.00 to -2.87
-3.00 to -3.87
-4.00 to -4.87
-5.00 to -5.87
-6.00 up

Lens Cylinder

reduce by .25D
reduce by .50D
reduce by .75D
reduce by 1.00D
reduce by 1.25D

Flat "K" of 43.00 & flatter -
Flat "K" of 43.00 to 46.00 -

use the 8.9mm Base Curve - 14.0 mm Diameter
use the 8.6mm Base Curve - 14.0 mm Diameter

NOTE:

Lens fitting should be conducted in a well illuminated room and overrefractions conducted using trial frame lenses.

Always Ahead of the Curve



800-223-1858

www.metro-optics.com