

MagniClear®

MagniClear®
plus

MagniClear® and MagniClear®*plus* are the first presbyopic lenses to combine the clear distance optics of a spherical lens, the intermediate visual performance of a multifocal and the predictably crisp, near power of a bifocal—all in one easy step.

MagniClear® provides:

- a wide range of posterior surface design options, including spherical, toric and eccentric

MagniClear®
plus provides:

- enhanced near power effect while controlling lens mass and peripheral edge lift

Parameters Available:

Diameters: 9.0, 9.2, 9.5, 9.6, & 10.0

Base Curves: Virtually any base curve your patient requires in .05mm steps

Distance Powers: +20.00 to -20.00 in 0.25 D steps

Add Powers: Up to +4.50 in 0.25 D steps

Our revolutionary front-surface design process has created a **new generation of presbyopic correction** that allows your patients to see **far, intermediate and near!**

Featured in

Bausch & Lomb
Boston®
MATERIALS

ARTOptical
contact lens, inc.

The Healthy Choice for Crisp Vision

toll-free ordering **800.253.9364**
consultation direct **800.566.8001**

toll-free fax **800.648.2272**

online **www.artoptical.com**

MagniClear® Plus

Fitting in 4 Simple Steps:

MagniClear®

MagniClear®
plus

FITTING & REFERENCE GUIDE

1. Base Curve Selection

Determine the base curve fit by referencing the flat K and the amount of corneal cylinder present.

<u>If Corneal Cylinder is:</u>	<u>Select Base Curve:</u>
Spherical to 1.25D.....	0.50 D steeper than K
1.50 to 2.25D.....	0.75D Steeper than K
2.50 to 3.00D.....	1.00D Steeper than K

2. Diameter Selection

Determine the diameter by the base curve selected.

<u>If Base Curve is:</u>	<u>Diameter Selection:</u>
8.50 to 8.45mm.....	10.0mm
8.40 to 8.20mm.....	9.6mm
8.15 to 7.50mm.....	9.5mm
7.45 to 7.20mm.....	9.2mm
7.15 to 6.90mm.....	9.0mm

3. Distance Power Selection

Determine the distance power by compensating for any vertex change (sphere powers of +/- 4.00 diopters or greater) and adjust for any tear layer change generated from going steeper than flat K.

4. Add Power Selection

Determine the near add power by adding 0.25D to the spectacle add power. If spectacle add power is +2.00, then contact lens near add power will be +2.25.

MagniClear/MagniClearplus Warranty:

Six month extended breakage and wettability warranty provides two refits per eye along with a worry free approach to fitting GP lenses, extended warranty to cover patient handling.

Featured in

Bausch & Lomb
Boston®
MATERIALS

MagniClear® and MagniClear® plus are
manufactured exclusively by:

ARTOptical
contact lens, inc.

toll-free ordering **800.253.9364**
consultation direct **800.566.8001**