Unique-pH®

Multi-Purpose Solution

Please read carefully and keep this package insert for future use in case you have a problem. **Unique-pH®** Multi-Purpose Solution cleans, conditions, disinfects and enhances initial lens wettability to promote comfortable wear of rigid gas permeable contact lenses. **Unique-pH®** Multi-Purpose Solution is a gentle, effective, chlorhexidine-free, thimerosal-free and benzalkonium chloride-free formulation for rigid gas permeable lens wearers. Always follow your eye care professional's instruction. Your eyes are unique. Your eye care professional may recommend additional products (such as OPTI-FREE® SUPRACLENS® Daily Protein Remover) or procedures to care for your lenses based upon your individual tear chemistry and lens wearing schedule. If your eye care professional recommends use of OPTI-FREE® SUPRACLENS® Daily Protein Remover, follow the directions for use on the OPTI-FREE® SUPRACLENS® package.

DESCRIPTION:

Unique-pH® Multi-Purpose Solution is a proprietary, sterile, aqueous solution buffered to approximate the same pH and tonicity of your eye. It contains hydroxypropyl guar, a proprietary wetting / conditioning polymer system, polyethylene glycol, Tetronic* 1304, boric acid, propylene glycol, and is preserved with POLYQUAD® (polyquaternium-1) 0.0011%, and edetate disodium 0.01%.

ACTIONS:

Cleans: **Unique-pH®** Multi-Purpose Solution removes dirt, protein deposits and debris from your contact lenses. No separate daily cleaner is required.

Conditions: This solution adjusts to your eye's natural tear pH to enhance the wettability and comfort of rigid gas permeable contact lenses upon insertion. This formula is designed to provide a soothing effect during lens wear.

Disinfects: This solution contains antimicrobial agents, which destroy harmful microorganisms (germs) commonly found on the surfaces of lenses.

Unique-pH® Multi-Purpose Solution may also be used as a diluent for OPTI-FREE® SUPRACLENS® Daily Protein Remover for simultaneous enzymatic cleaning during conditioning.

INDICATIONS (Uses):

Unique-pH® Multi-Purpose Solution is indicated for the cleaning, rinsing, disinfection, and conditioning of fluorosilicone acrylate and silicone acrylate rigid gas permeable contact lenses.

CONTRAINDICATIONS

(Reasons Not to Use):

If you are allergic to any ingredient in **Unique-pH®** Multi-Purpose Solution, do not use this product. **Unique-pH®** Multi-Purpose Solution is not for use with soft contact lenses.

WARNINGS:

PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYE. It is essential that you follow your eye care practitioner's directions and all labeling instructions for proper use of your lenses and lens care products including the lens storage case. EYE PROBLEMS, INCLUDING CORNEAL ULCERS CAN DEVELOP RAPIDLY AND LEAD TO LOSS OF VISION. DAILY WEAR LENSES ARE NOT INDICATED FOR OVERNIGHT WEAR AND SHOULD NOT BE WORN WHILE SLEEPING. IF YOU EXPERIENCE EYE DISCOMFORT, EXCESSIVE TEARING, VISION CHANGES, OR REDNESS OF THE EYE, IMMEDIATELY REMOVE YOUR LENSES AND PROMPTLY CONTACT YOUR EYE CARE PRACTITIONER.

It is recommended that contact lens wearers see their eye care practitioner twice each year or, if directed, more frequently.

- Studies have shown that smokers have a higher incidence of adverse reactions.
- Never touch the dropper tip of the container to any surface since this may contaminate the solution. Replace cap after use.
- To avoid contaminating your solution, do not transfer to other bottles or containers.

PRECAUTIONS:

- Never reuse this solution.
- Store at room temperature.
- Keep bottle tightly closed when not in use.
- Always wash and rinse your hands before handling your lenses.
- Your lenses must be cleaned and disinfected each time they are removed.
- Use before the expiration date marked on the container and carton.

- Keep this and all medications out of the reach of children.

ADVERSE REACTIONS (Problems and What to Do):
The following problems may occur:

- Eyes stinging, burning or itching (irritation)
- Excessive watering (tearing) of the eyes
- Unusual eye secretions
- Redness of the eyes
- Reduced sharpness of vision (visual acuity)
- Blurred vision
- Rainbows or halos around objects
- Sensitivity to light (photophobia)
- Dry eyes

If you notice any of the above problems, immediately remove and examine your lenses. If the problem stops and the lenses appear to be undamaged, thoroughly clean, rinse and disinfect the lenses and reinsert them. **Do not reinsert a damaged lens.** If the problem continues or a lens appears to be damaged, **IMMEDIATELY** remove your lenses and **IMMEDIATELY** consult your eye care practitioner.

If any of the above problems occur, a serious condition such as infection, corneal ulcer, neovascularization, or iritis may be present. Seek immediate professional identification of the problem and treatment to avoid serious eye damage. See your **Instructions for Wearers** booklet for more information.

GENERAL DIRECTIONS:

- Always wash and rinse your hands before handling your contact lenses.
- Always handle one lens, left or right, first in order to avoid mix-ups.
- Use only the solutions recommended by your eye care practitioner.
- Remove imprinted seal around bottleneck completely before use.

STEPS FOR PROPER DISINFECTION:

In order to provide appropriate care for your lenses, the following steps must be completed each time you remove your lenses. Do not skip any step.

Step 1: Cleaning

Rub both sides of the lenses carefully with 2 to 4 drops of **Unique-pH®** Multi-Purpose Solution. Be sure to clean each side of the lens.

Step 2: Rinsing

Thoroughly rinse with the **Unique-pH®** Multi-Purpose Solution or another rinsing solution recommended by your eye care practitioner.

Step 3: Conditioning and Storage

Place each lens in the correct side of the lens storage case, then fill each chamber with fresh **Unique-pH®** Multi-Purpose Solution completely covering the lenses. Close the case tightly. Store the lenses in **Unique-pH®** Multi-Purpose Solution for at least four hours prior to wearing your lenses. **Store your lenses in the unopened case until ready to wear, up to a maximum of 30 days.**

Before Wearing:

- Lenses may be removed from the storage case and placed directly into your eyes. For extra cushioning, a drop of the **Unique-pH®** Multi-Purpose Solution may be added to the inside of each lens prior to insertion.
- To prevent contamination and help avoid serious eye injury always empty and rinse your lens storage case with fresh solution and allow it to air-dry. Refill the case with fresh **Unique-pH®** Multi-Purpose Solution immediately before the next use of the case.

How Supplied:

Unique-pH® Multi-Purpose Solution is supplied in a sterile 3 fl. oz. (90 mL) size. The container and carton are marked with lot number and expiration date.

* Tetronic is a registered trademark of BASF Corporation.

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©2002 Alcon Laboratories, Inc.
*Tetronic es una marca registrada de BASF Corporation.
© 2002 Alcon Laboratories, Inc.