

IMPORTANT INFORMATION! PLEASE READ ALL OF THE FOLLOWING INFORMATION BEFORE USING CONTACT LENSES.

Contact Lenses shall be worn only with the assistance of qualified eye care practitioner's supervision.

Do not wear these lenses if following conditions are present:

- Acute and sub-acute inflammation or infection of the anterior chamber of the eye.
- Any eye disease, injury, or abnormality that affects the cornea, conjunctiva, or eyelids.
- Severe insufficiency of lacrimal secretion (dry eyes).
- Corneal hypoesthesia (reduced corneal sensitivity), if not-shakic.
- Any systemic disease that may affect the eye or be exaggerated by wearing contact lens.
- Allergic reactions of ocular surfaces or adnexa that may be induced or exaggerated by wearing contact lens or use of contact lens solutions.
- Any active corneal infection (bacterial, fungi, or viral).
- If eyes become red or irritated.
- Patients unable to follow lens care regimen or unable to obtain assistance to do so.

Contact lenses shall be prescribed by the eye care practitioner.

WARNING!

PROBLEMS WITH CONTACT LENS AND LENS CARE PRODUCT COULD RESULT IN SERIOUS INJURY TO THE EYE. IT IS ESSENTIAL THAT YOU FOLLOW YOUR EYECARE PRACTITIONER'S DIRECTION AND ALL LABELING INSTRUCTIONS FOR PROPER USE OF LENS AND LENS CARE PRODUCTS, INCLUDING THE LENS CASE.

EYE PROBLEMS, INCLUDING CORNEAL ULCERS, CAN DEVELOP RAPIDLY AND LEAD TO LOSS OF VISIONS; THEREFORE, IF YOU EXPERIENCE EYE DISCOMFORT, EXCESSIVE TEARING, VISION CHANGE, OR REDNESS OF THE EYE, IMMEDIATELY REMOVE YOUR LENS AND PROMPTLY CONTACT YOUR EYECARE PRACTITIONER.

IMPORTANT!

The contact lenses must be cared for properly. Please remember to follow the following guidelines before wearing the contact lenses.

1. Always wash your hands before inserting, removing or handling your contact lenses.
2. Use only the lens disinfection regimen recommended by your eye care practitioner.
3. Never use tap water.
4. Never share your lenses with anyone. This can lead to infection and the spread of disease.
5. Never sleep with lenses on. Wear only the number of hours recommended by your eye care practitioner.
6. Have your eyes examined yearly.
7. Do not let the hair spray or other cosmetic aerosols come in contact with your lenses.
8. Do not apply eyewash while wearing your lens.
9. Insert your lenses before putting makeup on.
10. If you experience any unusual redness, blurriness, discomfort, excessive itching or pain, remove your lenses immediately and consult your eye care practitioner.
11. Remember regular checkups and always follow the advice of your eye care practitioner.

INSERTION:

Remove the lens from the package. Make it a habit of handling the right lens first, which will prevent you from mixing up your lenses. Place the lens on your index finger. Check that the lens is clean and the right side up. Hold your upper lid open with your left hand's index and middle finger. At the same time, hold your lower lid open with the middle finger of your right hand. Looking straight at your eye in the mirror, slowly and gently place the lens on your eye. Then carefully release your eyelids and blink. You can gently massage your eyelid to center the lens.

REMOVAL:

Look up and hold lower lid down. Pull the lens gently downward with your index finger and thumb. Pinch the edge of the lens and remove it from your eye.

DISINFECTION:

Before inserting the lens after removing it from the eye, always rinse the lens with the STERILE contact lens solution suggested by your eye care practitioner. After removing your lenses, store them in a case filled with fresh STERILE contact lens solution suggested by your eye care practitioner.

These contact lenses are 45% Methafilcon A and 55% water and may be kept for up to 1 month after opening the "blister" package containing the lens and stored in proper contact lens case with STERILE contact lens solution suggested by your eye care practitioner.