

# Essential Xtra Fitting Guide

Lenses>SpecialtyGP>Multifocal> ESSentialXtraGP>Fitting Guide

## Parameter Selection:

### 1. Initial Base Curve/Diameter Selection

Select initial base curve according to the base curve selection chart.

Select diameter according to diameter selection chart.

## Base Curve Selection Chart

*Select base curve according to corneal cylinder*

Determine Flat K	0.00D to 0.62D	0.75D and UP
40.00 to 40.37	8.30	8.20
40.50 to 40.87	8.20	8.10
41.00 to 41.37	8.10	8.00
41.50 to 41.87	8.00	7.90
42.00 to 42.37	7.90	7.80
42.50 to 43.00	7.80	7.70
43.12 to 43.62	7.70	7.60
43.75 to 44.25	7.60	7.50
44.37 to 44.87	7.50	7.40

45.00 to 45.50	7.40	7.30
45.62 to 46.12	7.30	7.20
46.25 to 46.75	7.20	7.10
46.87 to 47.37	7.10	7.00

### Diameter selection chart

Select diameter according to corneal cylinder

Base curve	Myopes	Hyperopes
$\leq 7.30\text{mm}$	9.2mm	9.5mm
7.40 to 7.90mm	9.5mm	9.5mm
$\geq 8.00\text{mm}$	9.5mm	9.7mm

## 2. Add Series Selection

Select add series according to patient's add

Patient's Add	Add series
+0.50 to +1.50	Series 1
+1.75 to +2.25	Series 2

### 3. Lens Position and Movement

Evaluate lens position and movement; the ideal fit will be superior central (upper lid attachment) with a fluorescein pattern that demonstrates alignment along the flattest corneal meridian. Make base curve and diameter changes accordingly (see Troubleshooting guide).

### 4. Lens Power

Perform your over-refraction with loose trial lenses to determine the final distance Rx. Expect final Rx to be -0.50 D more than the existing contact lens Rx. Place the over-refraction in a trial frame and evaluate the transition from distance to near vision. If the over-refraction leads to acceptable distance but unacceptable near vision, reassess your base curve and/or add selection.

***Important Note:*** In order to maximize the ADD available in each series, the lens needs to translate upward along the vertical corneal meridian as the patient looks from distance to reading tasks. An upper lid attachment will facilitate the upward transition of the lens.

## Troubleshooting Guide

Clinical Observations	Common Symptoms	Fitting Goal	Remedy
<b>High-Riding Lens</b>	Excellent near,  Excessive minus  Over-refraction to achieve good  distance.	Lower the lens to pick up more distance power in the center of the lens.	Steepen B.C. by 0.10mm and specify a thin edge design  <b>And/Or</b>  Reduce lens diameter by 0.3mm, and specify a thin edge design  <b>Last Alternative</b> Steepen B.C. by 0.10, reduce diameter by

0.3mm, specify a thin edge design and increase the lens center thickness by 0.05mm

### **Lens Rides Temporally**

Near is usually excellent, distance requires additional -0.50D to -0.75D.

Get lens to center more by changing the dynamics between the upper lid and the lens.

**Diameter  $\leq$  9.3mm**

Increase lens diameter by 0.50mm and specify a minus carrier.

**Diameter > 9.3mm**

Increase lens diameter by 0.3mm and specify a minus carrier

### **Low-riding lens**

Good to excellent distance, near almost non-existent.

Need to achieve an upper lid attachment, making the greater plus power on the periphery of the lens available for near tasks.

**Diameter  $\leq$  9.3mm**

Increase lens diameter by 0.5mm and specify a minus carrier

**Diameter > 9.3mm**

Increase lens diameter by 0.3mm and specify a minus carrier

**Last Alternative**

Flatten B.C. by 0.10mm (only if the fit allows to do so).

<b>Lens that centers</b>	Excellent distance, reading difficult	Get the lens to translate upward upon downward gaze.	Increase lens diameter by 0.3mm and specify a minus carrier.
<b>Lens Rides Nasally</b>	Distance is difficult and near vision is almost non-existent.	Get the lens to center more, by changing the dynamics between the upper lid and the lens.	<p><b>Diameter <math>\leq</math> 9.3mm</b></p> <p>Increase lens diameter by 0.50mm and specify a minus carrier</p> <p><b>Diameter <math>&gt;</math> 9.3mm</b></p> <p>Increase lens diameter by 0.30mm and specify a minus carrier.</p>