clariti® 1day toric



Silicone hydrogel, daily disposable, toric contact lenses.



Features	Benefits	
High oxygen transmissibility.	Allows more oxygen to the cornea compared to hydrogel lenses, helping promote ocular health.	
Non surface treatment technology.	Patented AquaGen $^{\text{IM}}$ process ensures a highly wettable and lubricious lens surface promoting all-day comfort.	
Low modulus.	Soft and flexible material can help adaption from conventional hydrogel lenses.	
High water content.	Combined with AquaGen technology supports excellent all-day comfort.	
UVA and UVB protection.	Helps protect the eye from harmful UV rays.	
Aspheric back surface and blended edge design.	Optimizes lens movement on the eye. Minimises lens/conjunctival interactions.	
Aspheric Optics to minimise aberration.	Consistent and stable visual acuity.	
Smooth transitions between zones.	Superior comfort and handling for the patient.	
World's largest parameter range.	Fit a wider range of patients.	
High levels of refitting success.	Optimized design ensures interchangeability with other disposable lenses.	
Orientation mark.	Single orientation mark at 6 o'clock to allow simple assessment of fit.	



Product Specifications:

Water content 56% Base curve 8.60mm Diameter 14.3mm Centre thickness (@ -3.00DS) 0.105mm Dk/t (@ -3.00DS) 57 Modulus 0.5MPa Packaging 5 pack and 30 pack	Material	somofilcon A
Diameter 14.3mm Centre thickness (@ -3.00DS) 0.105mm Dk/t (@ -3.00DS) 57 Modulus 0.5MPa	Water content	56%
Centre thickness (@ -3.00DS) 0.105mm Dk/t (@ -3.00DS) 57 Modulus 0.5MPa	Base curve	8.60mm
Dk/t (@ -3.00DS) 57 Modulus 0.5MPa	Diameter	14.3mm
Modulus 0.5MPa	Centre thickness (@ -3.00DS)	0.105mm
	Dk/t (@ -3.00DS)	57
Packaging 5 pack and 30 pack	Modulus	0.5MPa
	Packaging	5 pack and 30 pack

Power Range

SPHERE POWER		CYLINDER POWER	AXIS
Plano to -6.00 -6.50 to -9.00	(0.25DS steps) (0.50DS steps)	-0.75 -1.25 -1.75	10° 20° 60° 70° 80° 90° 100° 110° 120° 160° 170°
+0.25 to +4.00	(0.25DS steps)		20° 70° 90° 110° 160° 180°
Plano to -6.00 -6.50 to -9.00	(0.25DS steps) (0.50DS steps)	-2.25	10° 20° 90° 160° 170° 180°