

## Bi-Toric Design From "K" & Rx

Example: K = 41.00 / 45.00 x 85 Rx = +0.50 -4.75 x 175

### » Determine the amount of corneal cylinder.

Answer: 4.00D Corneal Cylinder

### » Use the following chart for your base curve:

Corneal Cylinder	Base Curve	
	(Flat Meridian)	(Steep Meridian)
1.50 to 1.87D	Fit on K	Fit on K
2.00 to 2.37D	Fit on K	0.25D Flatter
2.50 to 2.62D	Fit on K	0.50D Flatter
2.75 to 5.87D	0.25D Steeper	0.50D Flatter
6.00D or greater	0.50D Steeper	1.00D Flatter

Answer: We will fit our base curve 0.25 steeper than the flattest meridian and 0.50 flatter than the steepest meridian.

Base Curve: 41.25 / 44.50

### » Put the spectacle Rx on the optical cross.



Vertex both meridians (using 12mm vertex distance).

Flat Meridian  $+0.50 \times 175 = +0.50$

Steep Meridian  $-4.25 \times 085 = -4.00$

### » Using SAMFAP, determine the power in both meridians. The flat meridian usually has the most plus or least minus in the prescription.

Flat Meridian +0.50

Steep Meridian -4.00

Steeptened Flat Meridian By -0.25

Flattened Steep Meridian By +0.50

Power in Flat Meridian +0.25

Power in Steep Meridian -3.50

Meridional Power = +0.25 / -3.50

Cylinder Power = +0.25 -3.75

### » Choose the diameter the same way you would for a spherical design using the average of the two base curves.

Base Curve	Recommended Diameter
39.00 to 43.87	9.6
44.00 to 45.87	9.3
46.00 to 48.37	9.0
48.50 to 50.00	8.7

**FINAL PRESCRIPTION:**

Base Curve 9.6:

41.25/44.50

Power:

+0.25 -3.75

Diameter:

9.6